

# Comprehensive Tumour Profiling

## CH/EEA/EU




Caris Life Sciences' comprehensive molecular profiling approach to assess DNA and RNA reveals a molecular blueprint to guide more precise and individualized treatment decisions.



### DNA

#### Whole Exome Sequencing

*SNVs, Indels, Copy Number Alterations, Karyotyping*



### RNA

#### Whole Transcriptome Sequencing

*Fusions, Variant Transcripts and Gene Expression*

## Technical Specifications

**Sufficient tumour content (>20% tumour nuclei) must be present to complete all analyses.**

Technical Information	MI Exome™ (Whole Exome Sequencing – DNA)	MI Transcriptome™ (Whole Transcriptome Sequencing – RNA)
Sample Requirements	FFPE block or 10 unstained slides with a minimum of 20% malignant origin for DNA and 10% malignant origin for RNA. Needle biopsy is also acceptable (4-6 cores).	
Tumour Enrichment (when necessary)	Microdissection to isolate and increase the number of cancer cells to improve test performance and increase the chance for successful testing from small tumour samples	
Number of Genes	22,000 genes	22,000 genes
Average Depth of Coverage (DNA) Average Read Count (RNA)	1,000x for 719+ clinical and research genes and 400-500x for all other genes	60 million
Positive Percent Agreement (PPA)	>95% for base substitutions at ≥ 5% mutant allele frequency; >99% for indels at ≥ 5% mutant allele frequency; >95% for copy number alterations (amplifications ≥ 6 copies)	>97%
Negative Percent Agreement (NPA)	>99%	>99%
Viruses	HPV 16 & 18 (Head & Neck, Anal, Genital, CUP) EBV (Head & Neck, Esophagogastric Junction, Gastric Adenocarcinoma, CUP. If positive result in cancer type not listed here, EBER ISH reflex to confirm EBV result.) MCPyV (Merkel Cell, Neuroendocrine – Poorly Differentiated (High-Grade)/Large or Small Cell Carcinoma, CUP)	
Genomic Signatures/Other	Loss of Heterozygosity (LOH) Microsatellite Instability (MSI) Tumour Mutational Burden (TMB) Human Leukocyte Antigen (HLA) Genotype Caris FOLFIRStai™ Caris GPSai™	


# Caris Molecular Testing – Complete Gene Coverage

As the pioneer in precision medicine, Caris was the first to provide WES and WTS for every patient. All molecular profiling orders include next-generation sequencing of all 22,000 genes.



## Whole Exome Sequencing (WES) DNA

- 22,000 genes
- 1,000x for 719+ clinical and research genes and 500x for all other genes
- SNVs, indels, copy number alterations and karyotyping
- 250,000 evenly-spaced genomic SNP
- Viruses
- Genomic signatures:
  - Loss of Heterozygosity (LOH)
  - Microsatellite Instability (MSI)
  - Tumor Mutational Burden (TMB)
- Other:
  - HLA Genotype



## Whole Transcriptome Sequencing (WTS) RNA

- 22,000 genes
- 60 million read count
- Gene fusions, variant transcripts and gene expression
- Novel translocation detection independent of intronic breakpoint

## Gene List

Listed below are the genes most commonly associated with cancer. Full 22,000 genes search is available on CarisLifeSciences.com.

ABL1	BCR	CSF1R	FANCC	GNAQ	LZTR1	MUTYH	PIK3CB	RAD51D	SOCS1
ABL	BLM	CTNNA1	FANCD2	GNAS	MAML2	MYB	PIK3R1	RAD54L	SPEN
ACVR1	BMPRI1A	CTNNB1	FANCE	H3F3A	MAP2K1	MYC	PIK3R2	RAD54L	SPPO
AIP	BRAF	CXCR4	FANCF	H3F3B	MAP2K2	MYCN	PIM1	RASA1	SRC
AKT1	BRCA1	CYLD	FANCG	HDAC1	MAP2K4	MYD88	PKN1	RB1	SSBP1
AKT2	BRCA2	CYP17A1	FANCI	HIST1H3B	MAP3K1	NBN	PMS1	RELA	STAG2
AKT3	BRD3	DDR2	FANCL	HIST1H3C	MAPK1	NF1	PMS2	RET	STAT3
ALK	BRD4	DICER1	FANCM	HNF1A	MAPK3	NF2	POLD1	RHOA	STK11
AMER1	BRIP1	DNMT3A	FAS	HOXB13	MAST1	NFE2L2	POLD2	RNF43	SUFU
APC	BTX	EGFR	FAT1	HRAS	MAST2	NFKBIA	POLD3	ROS1	TERT
AR	CALR	EGFRV8	FBXW7	IDH1	MAX	NOTCH1	POLD4	RPA1	TET2
ARAF	CARD11	EGLN1	FGFR1	IDH2	MED12	NOTCH2	POLE	RPA2	TFE3
ARHGAP26	CASP8	ELF3	FGFR2	INSR	MEF2B	NPM1	POLQ	RPA3	TFEB
ARHGAP35	CBFB	EP300	FGFR3	IRF4	MEN1	NRAS	POT1	RPA4	THADA
ARID1A	CCND1	EPHA2	FGFR4	JAK1	MET	NRG1	PPARG	RSPQ2	TMEM127
ARID2	CCND2	ERBB2	FGR	JAK2	MET Exon 14 Skipping	NSD1	PPP2R1A	RSPQ3	TMPRSS2
AR-V7	CCND3	ERBB3	FH1	JAK3	MGA	NTHL1	PPP2R2A	RUNX1	TNFAIP3
ASXL1	CD274	ERBB4	FLCN	KDM5C	MGMT	NTRK1	PRDM1	SDHA	TNFRSF14
ATM	CD79B	ERCC2	FLT1	KDM6A	MITF	NTRK2	PRKACA	SDHAF2	TP53
ATR	CDC73	ERG	FLT3	KDR	MLH1	NTRK3	PRKAR1A	SDHB	TRAF7
ATRX	CDH1	ESR1	FLT4	KEAP1	MLH3	NUMBL	PRKCA	SDHC	TSC1
AXIN1	CDK12	ETV1	FOXA1	KIF1B	MPL	NUTM1	PRKCB	SDHD	TSC2
AXIN2	CDK4	ETV4	FOXL2	KIT	MRE11	PALB2	PTCH1	SETD2	U2AF
AXL	CDK6	ETV5	FUBP1	KLF4	MSH2	PARP1	PTEN	SF3B1	VHL
B2M	CDKN1B	ETV6	FYN	KMT2A	MSH3	PBRM1	PTPN11	SMAD2	WRN
BAP1	CDKN2A	EWSR1	GALNT12	KMT2C	MSH6	PCNA	RABL3	SMAD4	WT1
BARD1	CHEK1	EXO1	GATA3	KMT2D	MSMB	PDGFRA	RAC1	SMARCA4	XPO1
BCL2	CHEK2	EZH2	GLI2	KRAS	MST1R	PDGFRB	RAD50	SMARCB1	XRCC1
BCL9	CIC	FANCA	GNA11	LCK	MTOR	PHOX2B	RAD51B	SMARCE1	XRCC2
BCOR	CREBBP	FANCB	GNA13	LYN	MUSK	PIK3CA	RAD51C	SMO	YES1



To order or learn more, visit [www.CarisLifeSciences.com](http://www.CarisLifeSciences.com).

US: 888.979.8669 | [CustomerSupport@CarisLS.com](mailto:CustomerSupport@CarisLS.com)

Intl: 00 41 21 533 53 00 | [InternationalSupport@CarisLS.com](mailto:InternationalSupport@CarisLS.com)

Where Molecular Science Meets Artificial Intelligence.

©2023 Caris MPI, Inc. All rights reserved. TNEU0001 v2 22 NOV 2023

